ENGLISH	CURAGEL®
Importer: Curamed <small>ophthalmics</small> Curamed Ophthalmics BV Nedereindeweg 405, 3437 PL Nieuwegein The Netherlands E.mail: info@curamed.nl, Web: www.curamed.nl NL-IM-000023059	

STERILE VISCOELASTIC SOLUTION FOR INTRAOCULAR INJECTION

Packaging:
 hydroxypropylmethylcellulose (HPMC) intraocular solution. Sterile syringe closed by a plunger stopper and a chlorobutyl tip cap, packaged in a pouch bag and dispatched in a box with sterile cannula.

Formula Each 1 ml of CuraGel contains:

Hydroxy propyl methyl cellulose (HPMC) IP-2% w/v	Calcium Chloride 2H ₂ O IP	-0.048% w/v
Sodium Chloride IP	Magnesium Chloride 6H ₂ O IP-0.03% w/v	
Sodium Acetate IP	Acetic Acid 1% Solution qs	
Sodium Citrate IP	Citric Acid 0.1465% Solution qs	
Potassium Chloride IP	Purified Water IP qs	

Properties:
 CURAGEL® is a viscoelastic, transparent, isotonic, sterile and apyrogenic solution of hydroxypropylmethylcellulose used as an adjuvant in ocular surgery. Injected into the anterior chamber of the eye, it provides a mechanical protection of the sensitive tissues of the eye, more specifically of the corneal endothelium. CURAGEL® maintains the transparency of the various parts of the eye and the depth of the anterior chamber during surgery. It helps to the implantation of the intraocular lens. CURAGEL® is inert and non toxic for any living tissues. It does not provoke any inflammation nor does it have any immunogenic, teratogenic or mutagenic characteristics.

Elimination of the product:
 CURAGEL® must be thoroughly removed after surgery by rinsing with irrigating solution. However, traces that may remain after surgery will almost totally dissipate (98%) through Schlemm's canal within 24 hours.

Indications:
 Surgical intraocular surgery of the anterior segment.

Precautions:
 Precautions are limited to those normally associated with the surgical procedures being performed. Instilling excessive amounts of CURAGEL® into the anterior or posterior chamber of the eye may increase intraocular pressure (IOP) temporarily. An excessive IOP may also be caused by a pre-existing glaucoma condition, by compromised outflow and by operative procedures and sequelae. Because these factors vary from case to case and are difficult to predict, the following precautions are recommended. Do not overfill the eye with CURAGEL. At the end of surgery, all remaining CURAGEL® should be removed by irrigation and/or aspiration. Carefully monitor the IOP especially during the immediate post-operative period. If a significant increase in pressure is observed, treat appropriately. Also: Avoid the injection of an air bubble at the end of the surgery. Do not use if sterile pouch is damaged. Do not use after the expiration date. Single use only. Dosage: Dosage varies with type of surgery. Usually a dose of 0.5 to 1 ml is injected into the anterior segment.

Ophthalmic surgical application:
 Cataract surgery: before lens extraction CURAGEL® is introduced into the anterior chamber to protect the corneal endothelium and to maintain a deep anterior chamber. Additional amounts may also be introduced prior to the implantation of an intraocular lens, and to coat the artificial lens and the surgical instruments.

Glaucoma surgery:
 CURAGEL® is instilled before trabeculectomy in order to reconstitute the anterior chamber. At the end of surgery, further instillation of CURAGEL® may be required to achieve good subconjunctival filtration and to prevent tissue adhesences.

Corneal trasplant surgery:
 prior to cornea trepanation, the anterior chamber is filled with CURAGEL. The donor graft should be placed into the proper position without removing CURAGEL. Then suturing can begin. Additional amounts of CURAGEL® may also be used in the anterior chamber of the donor eye to protect the corneal endothelium cells of the graft prior to trepanation, and to protect the exposed endothelium layer of the donor button during preparation of the recipient eye.

Incompatibilities:
 None currently known.

Adverse reactions:
 Transient rises of post-operative IDP have been reported in some cases. It could be controlled using appropriate medication.

Interactions:
 None currently known.

Storage instructions:
Store between 2°C and 35°C
Shelf life:
 3 years.

Instructions for use:
 Remove the sterile pouch from the box. Peel off the pouch bag above the sterile tray. Drop the syringe on the sterile tray. Remove the tip cap and attach the sterile 23 G angled cannula. Make sure the cannula is tightly connected with the Luer-Lok® tip of the syringe. Apply gentle pressure to the plunger in order to repel air from the cannula and to verify the syringe is operating properly. The syringe is ready for use.

NOTE: during surgery it is necessary to release the pressure on the plunger BEFORE pulling the cannula out of the eye, in order to avoid a spring effect which may cause the introduction of air in the cannula and the syringe.

GERMAN	CURAGEL®
---------------	-----------------

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